


# Repair kit for diaphragm pump for 321 & 301



Item number T206-321-100  
contains pos. 1-14

| Pos: | Pcs. | Drawing number | Item                             |
|------|------|----------------|----------------------------------|
| 1    | 2    | T07-100-009    | Diaphragm, outer vition          |
| 2    | 2    | T07-100-010    | Diaphragm, inner teflon          |
| 3    | 1    | L70-604-2000   | Hose multi oil, 1/2" x 2000 mm   |
| 4    | 2    | L70-604-425    | Hose multi oil, 1/2" x 425 mm    |
| 5    | 1    | L70-604-180    | Hose multi oil, 1/2" x 120 mm    |
| 6    | 1    | L63-42-20-12   | Ball bearing, 6004 RS 42x20x12   |
| 7    | 2    | L63-15-35-11   | Ball bearing, 6202 2RSH          |
| 8    | 1    | L92-115-000    | Locking ring outer, 15 mm        |
| 9    | 8    | L93-615-018    | O washer, 15,6 x 1,78 mm         |
| 10   | 10   | L99-118-12-02  | Rubber gasket, viton Ø13x19x2 mm |
| 11   | 4    | L99-513-000    | Ball Ø13, stainless              |
| 12   | 4    | T07-100-014    | Valve seat for pump 07           |
| 13   | 2    | L82-508-016    | Screw, 8 x16 mm                  |
| 14   | 1    | L85-15-21-01   | Washer, 15x21x1 mm               |

**BELKI Filtertechnik GmbH**  
Ostpromenade 54, D-52525 Heinsberg  
Tel +49 (0) 2452 10651-27  
www.belki.de • info@belki.de



Filtration with capacity

Hovedkontor / Headoffice / Hauptsitz:  
**BELKI teknik ApS**  
Fårupvej 19, DK-8840 Rødkærsbro  
Tel +45 8665 8811  
www.belki.dk • post@belki.dk

Presented by / Überreicht durch: